



OWNER'S MANUAL SHOULDER STRAP #140-SH

NOTE

- The #140-SH Shoulder Strap, supplements the standard #140 Patient Restraint Straps.
- Install only on Litters with the Shoulder Strap mounting crossbars as seen in Figure 1.
- The Red and Black straps of #140 are not utilized and may be discarded.
Important: Black strap from #140 is necessary if you are using flotation.




FIGURE 1

LIFESAVING SYSTEMS CORP.

220 ELSBERRY RD.
APOLLO BEACH, FL 33572
PHONE: 813/645-2748
FAX: 813/645-2768
www.lifesavingsystems.com