




**#406-1 Skid Shoe
Installation Instructions.**

Figure 1.

Place skids at head section, foot section, and middle of stretcher, as noted as black solid lines on figure 1. Clean skids with denatured alcohol. Remove backing from tape on shoe and firmly press shoe onto skid.

LIFESAVING SYSTEMS CORP.
220 ELSBERRY RD.
APOLLO BEACH, FL 33572
PHONE: 813/645-2748
FAX: 813/645-2768
www.lifesavingsystems.com
Lifesaving Systems Corp.
#F-750-406-1 O/M (Rev A) Updated 2014-07-14



LIFESAVING SYSTEMS CORP.
220 ELSBERRY RD.
APOLLO BEACH, FL 33572
PHONE: 813/645-2748
FAX: 813/645-2768
www.lifesavingsystems.com